# ReadyCal - Kit Anleitung

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# ReadyCal - Kit Manual

PSS ReadyCal mixtures are the easiest way to perform GPC/SEC calibrations. Use ReadyCals to obtain fast, precise and reproducible calibrations. The exact weight in the vial is known allowing the concentration to be altered by varying the amount of added eluent.

Advantages	Applicability
Fast	Conventional calibration
Reliable	Universal calibration
Easy and safe handling	Instrument validation
Free of contamination	Calibration of detectors

**Tip:** If you use an internal flow marker, please add it directly to the solvent used to prepare the ReadyCal solutions.

## How to use ReadyCals:

1. Select the appropriate vial and unscrew the cap.
2. Add the correct volume of solvent via syringe or pipette etc.

**Tip:** When viscometer or light scattering detectors are used weigh out the added amount of solvent precisely to calculate the concentration.

3. Screw on the vial cap and wait for 2 - 5 hours. Gently swirl the vial to homogenize.

**Tip:** Slightly shaking the vial during this time will speed up the dissolution process.

4. Inject each solution separately and use  $M_p$  values to construct a calibration curve.

**Storage** Place the tightly capped vials in the original box in a dry, dark, cool area; e.g. a refrigerator (4 °C).

**Expiration date** *Vials*  
Please see kit label and/or calibration curve delivered with any kit.

*Solutions*  
Use only freshly prepared solutions and use them up within two days. Solutions may only be stored in the refrigerator temporarily for a short time period. Do not freeze the standard solutions; freezing can destroy the sample.