

## Ensure WinGPC traceability and compliance (P/N 899-0033)

### Description

This software training course is designed for WinGPC Compliance Pack administrators, IT professionals responsible for CSV (computer system validation), and users responsible for validation. It offers lectures and practical PC sessions given by an experienced polymer chemist with extensive knowledge about chromatography and WinGPC UniChrom. Participants should be familiar with the basic WinGPC features and functions. Previous attendance of the WinGPC UniChrom software training course (P/N 899-0029) is recommended, but not required.

After successful participation each attendee should be able

- to administrate the WinGPC UniChrom Compliance Pack
- to identify and apply settings reflecting the company policy
- to control audit trails
- to set user level rights
- to perform meaningful IQ/OQ for GPC/SEC systems



### Program

<b>08.50</b>	<b>Login and Communication's Check</b>
<b>09.00</b>	<b>Welcome, speaker introduction, explanation of training tools</b>
<b>09.15</b>	<b>Introduction to WinGPC UniChrom with Compliance Pack</b> WinGPC UniChrom philosophy, Compliance Pack philosophy Guidelines for Macromolecular Chromatography Data Systems (MSDS) Overview national and international standards for GPC/SEC, electronic records, and others
<b>10.00</b>	<b>Installation and setup of the Compliance Pack</b> Installation requirements, infrastructure requirements User authentication (workstation/network options) Overview WinGPC Compliance Pack
<b>10.45</b>	<b>Compliance Pack configuration</b> The Compliance Pack Administration Tool User levels, user level rights, individual user rights Recommended configuration and settings
<b>11.15</b>	<b>Practical Session part I</b>
<b>12.00</b>	<b>Lunch</b>
<b>13.00</b>	<b>Efficient use of WinGPC UniChrom with Compliance Pack</b> User level dependent data acquisition and evaluation WinGPC projects with/without sample audit trail, option: reason for change Assign electronic signatures System audit trail, sample audit trail, instrument audit trail Traceability and compliance
<b>14.00</b>	<b>Practical Session part II</b>
<b>15.00</b>	<b>Comprehensive documentation</b> Use of the WinGPC UniChrom/PSS instruments Validation Binder PSS IQ/OQ tools, holistic test with the PSS EasyValid GPC/SEC validation kit Additional documentation available with the WinGPC ReportDesigner
<b>Approx. 16.00</b>	<b>Course review and end of training course</b>

